

Please click price below and you will be redirected to our pharmacy.

Pregnancy test strip



Pregnancy test strips - individually wrapped test devices. Each device contains an anti-alpha hCG capture antibody coated membrane and colloidal gold particles coated with mouse anti-beta hCG monoclonal antibody.

Active ingredient: [Pregnancy Test Strip](#)

Disease(s): Pregnancy / Fertility / Infertility

Pregnancy test strip 5mg

| | Per strip | Per Pack | Order |
|---------|-----------|----------|-----------------------------|
| 5 strip | \$ 0.60 | \$ 3.00 | Add to cart |

Accutest® Rapid Pregnancy Test

What is this product?

Accutest® Rapid Pregnancy Test - individually wrapped test devices. Each device contains an anti-alpha hCG capture antibody coated membrane and colloidal gold particles coated with mouse anti-beta hCG monoclonal antibody.

What should I know about Storage and Stability?

The test kit is to be stored refrigerated (2-8°C) or at room temperature (up to 30°C) in the sealed pouch for the duration of the shelf-life.

What should I know about Test Procedure?

Test device, patient's samples, and controls should be brought to room temperature (20-30°C) prior to testing. Do not open pouches until ready to perform the assay.

1. Remove the test device from its protective pouch (bring the device to room temperature before opening the pouch to avoid condensation of moisture on the membrane). Label the device with patient or control identification.

